

Alois

Maximum comfort minimum height

The Alois bed is designed for home or professional usage and is extremely stable due to its scissors lifting system. This nursing bed can be equipped with lowerable long wooden side rails or with metal side rails. The remote control, with soft-touch keys, takes care of all the adjustments.



VERMEIREN
we care for you

Comfort

The Alois is a low height bed for patients with Alzheimer's disease. Very low minimum height, maximum height up to 720 mm, and compatible with most patient lifts on the market.

Options

- 3 functions electrical bed with knee buckle
- Protective transport cover for the wooden panels
- Serum holder
- Remote control support
- Independent serum holder

Color

white



Technical data

■ Total width	1035 mm
■ Total height (without lifting pole)	660 - 1140 mm
■ Total length	2070 mm
■ Bedframe height (adjustable)	200 - 720 mm
■ Nominal load lifting pole	80 kg
■ Anti-Trendelenburg	no

Products designed to suit you

VERMEIREN always strives to come up with the best solution available to encompass your personal needs. As technology advances, so do we and as such we can offer you the most up to date products. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **Alois** complies with the highest norms, laid down by certified test institutes in most European countries.



Optional cradle for the remote

■ Adjustment of the back section	65°
■ Max. angle for lower leg	-16°
■ Max. thigh angle	30°
■ Weight	98 kg
■ Max. load (incl. patient, mattress and options)	170 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN GROUP

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com

For more information, contact your local dealer

