



# Goliat

## Heavy-duty rollator

The Goliat is a sturdy, foldable rollator that can carry a user up to 200 kg. The seat is filled with a foam cushion for optimal seating comfort. Furthermore, the user can adjust the height of the seat and the push handles to his or her preference. A standard hipbelt adds to the comfort and safety of the Goliat.

Total width	650 mm
Total length	750 mm
Seat height	490 - 560 mm
Seat width	470 mm
Seat depth	410 mm
Push handles	830 - 1020 mm
Multifunctional brakes with lock	
Weight	10,7 kg
Max. users weight	200 kg
Max. load basket	5 kg
Colour(s)	red/silver
Additional information	CE

Measurement tolerances:  $\pm 15$ mm / 1,5kg. We reserve the right to introduce technical changes. Models shown only for illustration purposes. Certain features may be optional.

