

KUARZ SYNC EON

Inspired by the past,
turning to the future



ABOUT US

VERMEIREN, a manufacturer since 1957, is a leading name in the medical equipment industry. Over the years, the Vermeiren Group has specialised in producing wheelchairs, nursing beds, shell chairs, comfort chairs, and other products designed to support patients and users in maintaining their independence and mobility.

Vermeiren's European development and design teams use state-of-the-art precision tools, ensuring reliability and quality in every product they create. Today, Vermeiren is a global group with over 1,300 employees and subsidiaries across Europe, the UK and Asia, supplying products to more than 45 countries worldwide.

Email : info@vermeiren.com
Website : <https://www.vermeiren.com>







MOBILITY

WITH DIGNITY AND STYLE

A modern electric scooter is more than just a way to get from A to B - it's freedom, dignity, and style rolled into one. With its sleek, modern design, it's the kind of ride you're proud to be seen on, whether you're meeting friends for coffee, heading out for a bit of shopping, or enjoying a leisurely afternoon in the park. It's about moving through life with confidence,

knowing you have the independence to go wherever you want, whenever you want.

This isn't just mobility - it's a lifestyle upgrade, combining comfort, practicality, and a touch of flair that turns every trip into a little celebration of you.





FAMILY

LINE UP

The new Vermeiren scooter collection offers a perfect balance of compact agility, refined versatility and premium comfort. The Eon, our smallest model, is designed for effortless indoor mobility and easy transport. The Sync bridges indoor and outdoor living with award-winning design, advanced comfort and everyday performance.

At the top of the range, the Kuarz delivers maximum stability, space and luxury for those who expect the very best. Three distinct personalities, one shared vision: redefining mobility with contemporary design, craftsmanship and freedom of movement.







EON

THE TIMELESS VOYAGER

Small in size, remarkable in impact. The Eon proves that true innovation doesn't need space to shine. Its ultra-compact design makes everyday mobility effortless whether navigating tight interiors or busy streets. Every curve serves a purpose, combining practicality with elegance. Recognized

for its award-winning design, the Eon delivers a seamless blend of convenience, comfort, and style in one refined, space-saving package.

The ideal scooter to increase mobility and bring back your freedom of movement.



SBS* LED FRONT LIGHTS



LED INDICATORS



FOLDING BACKREST



ADJUSTABLE ARMRESTS



LED REAR LIGHTS



ANTI-TIPPING



DRIVE OR NEUTRAL MODE



CONTROLLER WITH BATTERY INDICATION



STEERING MUDGUARDS

*SBS: See and Be Seen - Advanced LED light technology for increased safety when driving in low light conditions.



SYNC

THE HARMONIOUS CONNECTION

Bridging worlds with effortless elegance, the Sync is the perfect balance of versatility and design. Positioned between the ultra-compact Eon and the powerful Kuarz, it adapts seamlessly to both indoor comfort and outdoor freedom. Its award-winning design is not just about looks—it enhances

usability at every level. The Sync redefines what a mid-size scooter can be: intuitive, refined, and perfectly in tune with modern mobility.

You don't need to look any further as this scooter provides the best of both worlds



SBS* LED FRONT LIGHTS



LED INDICATORS



COMPACT STORAGE



ROTATING AND ADJUSTABLE SEAT



LED REAR LIGHTS



ANTI-TIPPING



DRIVE OR NEUTRAL MODE



CONTROLLER WITH BATTERY INDICATION



FOLDING BACKREST

*SBS: See and Be Seen - Advanced LED light technology for increased safety when driving in low light conditions.



KUARZ

DESIGNED WITH PRECISION

At the pinnacle of the range, the Kuarz combines commanding presence with exceptional comfort and performance. Designed for those who expect more from their mobility scooter, it delivers effortless stability, generous space and refined craftsmanship in every detail. Its contemporary, award-winning design is not only visually striking, but engineered to

enhance confidence and control in every journey. The Kuarz redefines premium mobility: powerful, sophisticated and uncompromisingly comfortable.

You don't need to compromise when this scooter delivers luxury, performance and freedom in perfect harmony.



SBS* LED FRONT LIGHTS



LED INDICATORS



COLLAPSIBLE CHAIR



PARKING BRAKE



LED REAR LIGHTS



ANTI-TIPPING



ROTATING AND ADJUSTABLE SEAT



LCD DISPLAY WITH ADAPTIVE BRIGHTNESS



PREMIUM SEAT FINISHING

*SBS: See and Be Seen - Advanced LED light technology for increased safety when driving in low light conditions.





EUROPEAN
PRODUCT
DESIGN
AWARD®
2025 WINNER



reddot winner 2026

All our 3 scooter models
Eon, Sync and Kuarz
have been awarded for its
outstanding design and new
approach towards scooters.





CUSTOMIZE YOUR RIDE

Options are what turn a scooter from “just transport” into your scooter. Mirrors, rearview cameras, and rain covers aren't just accessories—they're little touches that make every ride smoother, safer, and more enjoyable. A dust cover keeps it looking sharp, while a rollator holder means your mobility always comes with you.

Whether it's comfort, style, or convenience, customizing your scooter lets you match it perfectly to your daily life. After all, why settle for standard when you can have a ride that's as unique, practical, and ready-for-anything as you are?

From left to right: Quadri-Light rollator holder, MOVO rain/sun protection, set of mirrors, slick or studded tyres, colored hubcap matching your scooter color, 286i rollator holder, topcase at the back, shopping basket at the front, double crutch holder, storage protection cover, fully integrated rearview camera. Certain options may be tied to a specific scooter so please check our order forms for up to date information.






COLORS

OWN IT

All 3 models are available in a vibrant range of six colors, because your scooter should reflect you. Whether you want it to match to color of your eyes, coordinate with your wardrobe, or simply choose your all-time favorite shade, color is a powerful way to express your personality. A bold red might say "adventurous spirit," while a cool grey whispers "calm and confident."


From stylish neutrals to eye-catching brights, there's a perfect match for every mood and taste. After all, a scooter isn't just transportation, it's part of your personal style, and the right color makes all the difference.



Arctic Green




Burgundy Red



Cool Grey



Matt Gold



Olive Green



Pearl Orange



Cool Grey

Why choose Cool Grey for your Kuarz? Because it speaks of modern confidence, understated elegance, and a calm, dependable nature, ideal for those who prefer style that whispers rather than shouts.



Pearl Orange

Pearl Orange suggests a personality that's energetic, creative and full of optimism. It reflects a zest for life, a love of standing out, and an ability to inspire warmth and positivity in others.





Olive Green

Olive Green as a color for your Sync? This embodies harmony, grounded confidence and a connection to nature, ideal for those who appreciate calm strength and a style that feels both fresh and enduring.



Mat Gold

Mat Gold reflects a personality that values refinement, confidence and timeless elegance. You appreciate life's finer details, and enjoy making a statement that feels both subtle and distinguished.



Burgundy Red

Burgundy Red on the Sync scooter speaks of depth, warmth, and quiet strength. It's a choice for those who embrace life with passion, appreciate fine details, and exude an effortlessly classic charm.



Arctic Green

The Eon in Arctic Green reflects a calm, balanced personality with a love for nature and fresh beginnings - perfect for users who value harmony, optimism, and a refreshing, confident style.




TECHNIQUE

IN DEPTH

When it comes to scooters for seniors, the technical details aren't just numbers, they're the keys to freedom on four wheels. The width and length decide whether you glide gracefully through doorways or have to do a 2-point turn in the hallway. Weight tells you if it's a sturdy road companion or easier to nudge around.

Range means knowing if you can make it to the café and back without hunting for a charging plug. Even the top speed matters - fast enough to feel the breeze, slow enough to keep your hat on. These specifications make sure your ride is safe, comfy, and joyfully yours.



KUARZ

■ Total width	680 mm
■ Total length	1505 mm
■ Total height	1345 mm - 1395 mm
■ Folded height	870 mm
■ Seat height	480 mm - 530 mm
■ Seat width	525 mm
■ Seat depth	480 mm
■ Backrest height	550 mm
■ Footplate height	200 mm

■ Obstacle climbing	100 mm
■ Turning diameter	3400 mm
■ Batteries	80 Ah (C20), 12V DC, non-spillable AGM battery
■ Degree of protection	IPX4
■ Motors	750 W
■ Speed*	15 km/h
■ Driving range*	45 km
■ Weight	158 kg
■ Maximum user weight	150 kg

*The driving range is calculated according to European standards. *Speed and driving range may deviate (ca. 15%). Permitted dimensional and weight tolerances: 15 mm, 1.5 kg and 1.5°. Subject to technical alterations. Images for illustration purposes only. The batteries and speed maybe different in your country depending on local laws and regulations.

SYNC

■ Total width	630 mm
■ Total length	1235 mm
■ Total height	1275 mm - 1325 mm
■ Folded height	830 mm
■ Seat height	420 mm - 470 mm
■ Seat width	470 mm
■ Seat depth	495 mm
■ Backrest height	530 mm
■ Footplate height	185 mm

*The driving range is calculated according to European standards. *Speed and driving range may deviate (ca. 15%). Permitted dimensional and weight tolerances: 15 mm, 1.5 kg and 1.5°. Subject to technical alterations. Images for illustration purposes only. The batteries and speed maybe different in your country depending on local laws and regulations.

■ Obstacle climbing	110 mm
■ Turning diameter	3530 mm
■ Batteries 38 Ah or 50 Ah (C20), 12V DC, non-spillable AGM battery	
■ Degree of protection	IPX4
■ Motors	450 W / 500 W / 650 W
■ Speed*	6 km/h / 10 km/h / 12 km/h / 15 km/h
■ Driving range*	28 km - 43 km
■ Weight	109 kg
■ Maximum user weight	140 kg

EON

■ Total width	565 mm
■ Total length	1170 mm
■ Total height	1050 mm - 1100 mm
■ Folded height	900 mm
■ Seat height	460 mm - 520 mm
■ Seat width	470 mm
■ Seat depth	420 mm
■ Backrest height	480 mm
■ Footplate height	140 mm

*The driving range is calculated according to European standards. *Speed and driving range may deviate (ca. 15%). Permitted dimensional and weight tolerances: 15 mm, 1.5 kg and 1.5°. Subject to technical alterations. Images for illustration purposes only. The batteries and speed maybe different in your country depending on local laws and regulations.

■ Obstacle climbing	80 mm
■ Turning diameter	2700 mm
■ Batteries 38 Ah (C20), 12V DC, non-spillable AGM battery	
■ Degree of protection	IPX4
■ Motors	320 W
■ Speed*	6 km/h or 10 km/h
■ Driving range*	30 km
■ Weight	89 kg
■ Maximum user weight	140 kg



PAST MEETS THE FUTURE

In our two movies we strived to make a connection between the past and the future. Elegant lines and shapes are timeless, either in a classic mode of transportation that was built with love and care or in a new and trendy design for an electrical scooter. **Scan the QR codes on top of this page and be prepared to discover the feel and looks of our new scooters.**



SCAN AND
DISCOVER
THE KUARZ



SCAN AND DISCOVER
THE EON AND SYNC





VERMEIREN
we care for you



R.E.: **VERMEIREN** GROUP, Vermeirenpaleis 1/15, 2920 Kalmthout - Belgium
info@vermeiren.com - www.vermeiren.com - BL Kvarz Sync Eon World 2026-05 v1.0