



# V300 XL

## Maximised duties

The V300 XL is the top model in the V series.  
To reduce the overall weight and ensure  
easy transport, the V300 XL is made of  
lightweight aluminum.



## Heavy duty wheelchair

The **V300 XL** offers many technical adjustment options and can thus be individually adapted to the user. With a maximum user weight of 170 kg, the V300 XL is one of our lightest wheelchairs for people with a higher body weight as it only weighs 21 kg.

### Colours

C59

- ▶ 3 Seat widths
- ▶ 4 seat heights
- ▶ 2 seat depths
- ▶ 3 different armrest models
- ▶ maximum user weight: 170 kg
- ▶ Variety of different options

### Technical data

■ Total width	740 / 770 / 810 mm
■ Total height	840 - 930 mm
■ Total length	1085 mm
■ Folded width	350 mm
■ Seat width	530 / 560 / 600 mm
■ Height armrests	220 - 240 mm

## Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

### Quality guaranteed

Regarding quality, the **V300 XL** complies with the highest norms, laid down by certified test institutes in most European countries.



■ Seat height	440 - 530 mm
■ Height backrest	400 mm
■ Seat depth	490 - 530 mm
■ Weight	21 kg
■ Maximum user weight	170 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

## VERMEIREN Group N.V.

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

[www.vermeiren.com](http://www.vermeiren.com)  
[info@vermeiren.com](mailto:info@vermeiren.com)



For more information, contact your local dealer

**VERMEIREN**  
we care for you