

MAX. USER WEIGHT
250 KG

Forest 3+

bariatric wheelchair

The Forest 3+ is an electric wheelchair specially designed for large seat widths and for a higher user weight up to 250 kg. Complete with the latest electronic controller and a wide range of options.



Standard equipped with the LiNX controller.



FOREST 3+ in numbers

- Compact frame with a seat width from 500 mm
- Ground clearance up to 75 mm
- Stability axle ensuring contact of all four wheels with the ground on uneven terrains up to 100 mm
- Seat height 450 mm

V-STAB Technology influences the driving behavior of the power wheelchair to the extent that during driving maneuvers the direction of travel desired by the joystick can be mastered without joystick corrections regardless of uneven ground.

Colors (matte and glossy colors)



Options



SE09 - Additional controller



SE04 Chin controller



SE55 - Reinforced joystick



L58 - Multi-adjustable headrest



SE83 - Mini joystick



B66MAL Hemiplegic armrest

Technical data

■ Total width	705 / 720 / 760 / 820 / 870 mm
■ Total length	1145 / 1163 mm
■ Total height	1100 mm
■ Seat height	450 mm
■ Seat width	500 / 550 / 600 / 650 / 700 mm
■ Seat depth	450 - 550 mm

■ Batteries / Motors	2 x 12V / 80Ah - 350W
■ Speed	6 km/h* or 10 km/h*
■ Driving range	49 km*
■ Weight	134,35 kg
■ Maximum user weight	250 kg
■ Additional information	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer