

Eclips X4 Kids

Ultra light kids wheelchair

Constructed in aluminum, the Eclips X4 kids offers you an ultra light wheelchair combined with sturdiness and ease of use. Equipped with clothing protectors results in a very modern look. Without the wheels and the footrests, the Eclips X4 kids weighs no more than 9 kilograms, which translates in easy transport.



VERMEIREN
we care for you

Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.



Clothing protectors



Quickrelease rear wheels



Aluminum front forks adjustable in height



Internal and external collapsible footplates

Technical data

■ Total width	500 / 550 mm
■ Total height	820 - 980 mm
■ Total length	830 - 1000 mm
■ Folded width	310 mm
■ Seat width	320 / 370 mm
■ Height armrests	150 - 180 mm

Quality guaranteed

Regarding quality, the **Eclips X4 Kids** complies with the highest norms, laid down by certified test institutes in most European countries.

Colours



Model T30 push version with AR 12'' wheels

■ Seat height	420 - 470 mm
■ Height backrest	350 - 400 mm
■ Seat depth	310 - 370 / 330 - 390 mm
■ Weight	15 kg
■ Maximum user weight	80 kg
■ Additional info	☎ CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer

VERMEIREN
we care for you