

# 139E

## COMMODE CHAIR



total width	560 mm
seat depth	420 mm
seat width	420 mm
height armrest	226 mm
seat height	528 mm
total height	950 mm
total length	820 mm
backrest height	375 mm
weight	13 kg
max. load	120 kg
color frame	Blue
upholstery	Black

Deviations in sizes and weight (1,5 mm and 1,5 kg) are possible - Subject to technical alterations.



**N.V. Vermeiren N.V.** - Vermeirenplein 1-15 - B-2920 Kalmthout - Belgium

Tel.: 00 32 (0)3 620 20 20 - Fax: 00 32 (0)3 666 48 94 - info@vermeiren.be - visit our website: [www.vermeiren.com](http://www.vermeiren.com)