












Technical information

										
total width	seat width	height of armrests	seat height	total height	total length	seat depth	back height	angle back	users weight	weight
660	540	0-170	510	1200	750	440	780	0°- 45°	130	19 ^(P*) 20 ^(L*) 26 ^(N*)

* P: Provence, L: Languedoc, N: Normandie

Deviations in sizes and weight (15 mm and 1,5 kg) are possible



Seat without seams



Adjustable footrest



Normandie on wheels

Normandie adjustable back/legrest

Languedoc adjustable back

Provence fixed model



N.V. Vermeiren N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium


Phone: +32 (0)3 620 20 20
Fax: +32 (0)3 666 48 94

www.vermeiren.com
info@vermeiren.be



VERMEIREN
Opening up new horizons in mobility

For more information, contact your local dealer.



R.E.: N.V. Vermeiren N.V., Vermeirenplein 1/15 - 2920 Kalmthout - Belgium - 12/2010

Subject to technical alterations



Normandie



■ **Normandie** : adjustable back (45°) with a simultaneously moving legrest

■ **Languedoc** : adjustable back (45°)

■ **Provence** : fixed back

All models are equipped with armrests adjustable in height.


Adjustable back

Adjustable armrests

Adjustable legrest



Option: serumholder



Transport wheels



Option: adjustable height



Option: retractable footrest (version Languedoc and Provence)

Option: Electrical version

Provence

- Inclination of the back and legrests
- Long armrests for a complete support of the arm when the backrest is totally reclined

Option: Wide anatomical armrests

Languedoc



Skai Orange

Skai Light Green

Skai Light Blue

Skai Deep Green

Skai Dark Red

Skai Yellow

Skai Brown

Skai Night Blue

Tissue red

Tissue Blue