



\* Quantities are for one side  
 \* Only right side depicted

Material			
No.	Plan	Article	Quantity
1	X005982	2080000	3
2	X005983	2080100	1
3	M6x40 (DIN912/8.8)	1003909	3
4	M6x12	1008492	3
5	(Wheel 12")	1914618	1
6	M12x95 (DIN931/8.8)	3992074	1
7	XC00378	1C40254	1
8	EC00245	3991717 (L)	1
	EC00246	3991718 (R)	1
9	X001308	G051700	1
10	X001309	G051800	1

Left brake on the right frame side (3991717)  
 Right brake on the left frame side (3991718)

E001163 / C

<b>9300 SERIES T30 FIXED WHEEL</b>				Plan No.
				E001163 / C
Scale: 2:5	Drafted	Verified	Confirmed	
Date	26/10/2020	17/06/2021	17/06/2021	
Initials	TT	SD	PJ	