


Bill of Material			
No.	Plan	Article	Quantity
1	E00069095	2998116 (gray soft)	1
	E00069095	2998117 (black soft)	1
	E00069095	2998118 (cordura)	1
2	Y011148	1669754	6
3	Y011987	1669854	1
4	M8x35 with cam	1001526	3
5	Nut M8	1902024	3
6	M8 DIN126 Black	1996214	3
7	Cap 15x15	1005899	1
8	Star knob M8x16	2990024	1
10*	Y017183	1748154	1
11*	M6x12 ISO7380 Black	3992768	4
12***	EP00216	2998018 (D25.4)	1
		1913740 (D22)	1

\* use with L34 backrest

\*\* use with backrest pipes or with standard B90 pipe

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

				---
L55B				Plan No.
				E00069094/D1.1
Scale: 1:3	Drafted	Verified	Confirmed	
Date	18/10/2023	16/02/2024	16/02/2024	
Initials	PM	PM	AS	