


- * Left side only depicted
- * Other parts are kept from wheelchair
- * Quantities are for one pair

Material			
No.	Plan	Article	Quantity
1	Wheel 22" (PU)	1914625	2
	Wheel 22" (AIR)	1911073	2
2	Plastic fork	1992599	2
3	Wheel 140x37	1914613	2
4	M8	1998994	2
5	M12x1.25	1994605	2

E000636 / C

1907449 / 1908226
(PU)
(AIR)

S1 REPLACEMENT KIT				Plan No. E000636 / C
TYPE V200 / V300 / V500				
Scale: 1:6	Drafted	Verified	Confirmed	
Date	31/08/2011	09/12/2019	09/12/2019	
Initials	SD	TT	PJ	